

Mais confiança ao reforçar



Ressecção pulmonar



Ressecção gástrica



Gastrectomia em manga



Anastomose intestinal



Ligadura do apêndice atrial esquerdo

DESEMPENHO que vem da inovação

O Reforço GORE® SEAMGUARD® consiste em uma matriz exclusiva e patenteada de tecido sintético, cujo desempenho é constante e previsível. Sempre.

O produto foi concebido para reduzir a incidência de fistulas e sangramentos perioperatórios em diversos tipos de cirurgias abertas ou minimamente invasivas^{1,2}. O Reforço GORE® SEAMGUARD® é totalmente absorvido em seis a sete meses, eliminando o risco de uma resposta inflamatória prolongada.

- Maior quantidade de literatura clínica indicando menos fistulas
- Evidências mais fortes na literatura clínica indicando menos sangramentos
- Material sintético 100% bioabsorvível
- Material uniforme e fino
- Nenhum relato de migração nem de complicações de longo prazo
- Sem risco de contaminação de fontes animais

1. Franklin. Safety and efficacy of the use of bioabsorbable "Seamguard" in colorectal surgery at the Texas Endosurgery Institute. *Surgical Laparoscopy, Endoscopy & Percutaneous Techniques* 2005;15(1):9-13.

2. Consten ECJ, Gagner M, Pomp A, Inabnet WB. Decreased bleeding after laparoscopic sleeve gastrectomy with or without duodenal switch for morbid obesity using a stapled buttressed absorbable polymer membrane. *Obesity Surgery* 2004;14:1360-13 66.



SEAMGUARD®
BIOABSORBABLE

STAPLE LINE REINFORCEMENT



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)
00800.6334.4673 (Europe)
800.437.8181 (United States)
928.779.2771 (United States)

goremedical.com

Configurado para grampeadores endoscópicos

	GRAMPEADOR	COR DO CARTUCHO	NÚMERO DE CATÁLOGO
Endocirurgia Ethicon	ENDOPATH® ETS 45 para transecção endoscópica	Branco, azul, verde	1BSGETS45 / 12BSGETS45*
	ENDOPATH® ETS Flex 45	Branco, azul, verde	1BSGETS45 / 12BSGETS45*
	ENDOPATH® EZ 45	Azul, verde	1BSGEZ45 / 12BSGEZ45*
	ENDOPATH® ETS Flex 60	Branco, azul, verde	1BSGETS60 / 12BSGETS60*
	Grampeador Echelon 45 ENDOPATH®	Branco, azul, dourado, verde	1BSGEC45 / 12BSGEC45*
	Grampeador Echelon 60 ENDOPATH®	Branco, azul, dourado, verde	1BSGEC60 / 12BSGEC60*
	Echelon FLEX 60 ENDOPATH®	Branco, azul, dourado, verde	1BSGEC60A / 12BSGEC60A*
Covidien Autosuture	Endo GIA™ Universal 45 reto	Branco, azul	1BSGUSS45B / 12BSGUSS45B*
	Endo GIA™ Universal Roticulator 45	Branco, azul	1BSGUSS45B / 12BSGUSS45B*
	Endo GIA™ Universal 45 reto	Verde	1BSGUSS45G / 12BSGUSS45G*
	Endo GIA™ Universal Roticulator 45	Verde	1BSGUSS45G / 12BSGUSS45G*
	Endo GIA™ Universal 60 reto	Branco, azul	1BSGUSS60B / 12BSGUSS60B*
	Endo GIA™ Universal Roticulator 60	Branco, azul	1BSGUSS60B / 12BSGUSS60B*
	Endo GIA™ Universal 60 reto	Verde	1BSGUSS60G / 12BSGUSS60G*
	Endo GIA™ Universal Roticulator 60	Verde	1BSGUSS60G / 12BSGUSS60G*

* Doze pacotes por caixa

Configurado para grampeadores abertos

	GRAMPEADOR	NÚMERO DE CATÁLOGO
Endocirurgia Ethicon	Cortador linear PROXIMATE® 55 (Grampeadores: TLC55, TCT55; Recargas: TCR55, TRT55)	10BSGPROX55
	Cortador linear PROXIMATE® 75 (Grampeadores: TLC75, TCT75, TCD75; Recargas: TCR75, TRT75, TRD75)	10BSGPROX75
	Cortador linear PROXIMATE® 100 (Grampeadores: TLC10, TCT10; Recargas: TCR10, TRT10)	10BSGPROX100
Covidien Autosuture	GIA™ DST 60 (Grampeador: GIA6038S; Recarga: GIA6038L) (Grampeador: GIA6048S; Recarga: GIA6048L)	10BSGDST60
	GIA™ DST 80 (Grampeador: GIA8038S; Recarga: GIA8038L) (Grampeador: GIA8048S; Recarga: GIA8048L)	10BSGDST80
	GIA™ DST 100 (Grampeador: GIA10038S; Recarga: GIA10038L) (Grampeador: GIA10048S; Recarga: GIA10048L)	10BSGDST100

Configurado para grampeadores circulares

	GRAMPEADOR	NÚMERO DE CATÁLOGO
Endocirurgia Ethicon	ENDOPATH® 25 e Proximate 25 (Grampeadores: ECS25, CDH25, SDH25)	1BSGC25
	ENDOPATH® 29 e Proximate 29 (Grampeadores: ECS29, CDH29, SDH29)	1BSGC2829
	ENDOPATH® 33 e Proximate 33 (Grampeadores: ECS33, CDH33, SDH33)	1BSGC33
Covidien Autosuture	CEEA 25 (Grampeadores: 110214L, 111985)	1BSGC25
	CEEA 28 (Grampeadores: 110238L, 111987)	1BSGC2829
	CEEA 31 (Grampeadores: 110276L, 111989)	1BSGC31
	CEEA 34 (Grampeador: 111981)	1BSGC3334
	DST Series™ CEEA 25 (Grampeadores: EEA25, EEAXL25, EEA2535, EEAXL2535, EEAORVIL25)	1BSGC25
	DST Series™ CEEA 28 (Grampeadores: EEA28, EEAXL28, EEA2835, EEAXL2835)	1BSGC2829
	DST Series™ CEEA 31 (Grampeadores: EEA31, EEAXL31)	1BSGC31
	DST Series™ CEEA 34 (Grampeadores: EEA34, EEAXL34)	1BSGC3334

ETHICON, ENDOPATH, and PROXIMATE are trademarks of Johnson & Johnson Corporation.
AUTOSUTURE, ENDO GIA, and Tri-staple are trademarks of Covidien AG or its affiliates.

Products listed may not be available in all markets
GORE®, PERFORMANCE THROUGH INNOVATION, SEAMGUARD®, and designs are trademarks of W. L. Gore & Associates.
© 2011, 2012, 2013 W. L. Gore & Associates, Inc. AL2961-PT7 AUGUST 2013