



SOUND REACH[®]

Ultrasonic Surgical System

SOUND REACH provides speed, safety, ease of use and flexibility in surgeries where an advanced energy device is necessary.



Cincinnati, USA
4480 Lake Forest Dr, Suite 414
Cincinnati, OH 45242, USA
Tel: +1(513)769-1918 Fax: +1(513)769-1921

Tianjin, China
4th Floor, Bldg. B, No. 5 of Fourth
Street, TEDA, Tianjin, 300457, China
Tel: +86 22-25325428 Fax: +86 22-66200915

Beijing, China
Room 2015, 20 Floor, The Beijing Silver Tower,
No.2, East Third Ring North Road, Chaoyang
District, Beijing, 100027, China
Tel: +86 10-65675672 Fax: +86 10-65678258

<http://www.reachsurgical.com>

Please always refer to the Instructions for Use

May 16, 2018



Dissection, cutting and coagulation in a single device

Increases quickness and efficiency in surgery, minimizing the need to change instruments



No electricity flow through patients or medical practitioners

Provides safety during surgery



Minimal thermal spread

Reduces damage to adjacent tissue and organs



Reliable and secure vessel sealing across a wide range of vessels

Reduce the number of surgical devices needed to achieve hemostasis



One-handed activation and rotation

Makes it easy to use for surgeons



Same transducer for all shear sizes

Reduces costs



Shears compatible with other platforms

Provides savings and creates more value



* Sound Reach shears were shown to be able to securely and reliably seal 3-5mm vessels on porcine mesenteric vessels in clinical studies.

SOUND REACH SHEAR SIZES

CS1405H

- Curved 5mm
- 14cm long hand activated shear
- For open procedures, ENT, etc



CS2305H

- Curved 5mm
- 23cm long hand activated shear
- For pediatric, general surgeries



CS3605H

- Curved 5mm
- 36cm long hand activated shear
- For general, bariatric, gynecology, urology and colorectal endoscopic surgeries



OTHER INFORMAITON

DEVICES		
CODE	DESCRIPTION	MINIMUM ORDER QTY
CS1405H	14cm long shear	1 carton, 24pcs
CS2305H	23cm long shear	1 carton, 24pcs
CS3605H	36cm long shear	1 carton, 24pcs

CAPITAL		
CODE	DESCRIPTION	MINIMUM ORDER QTY
CSUS6000	CSUS6000 generator system (with transducer, foot switch and transducer cable)	1pc
TRA5-E	Transducer for CSUS6000 generator	1pc
TRA5 cable	Transducer cable	1pc
FSW1	Foot switch	1pc
CART6000	Generator cart	1pc